**HEMEDIX**  
English  
**QFlow 500 Probe™**  
For Use with the **Bowman Perfusion Monitor™**  
**INSTRUCTIONS FOR USE**

**1. Introduction**  
The QFlow 500 Probe incorporates a thermistor at the distal tip, which is controlled by the Bowman Perfusion Monitor, Model 500. The QFlow Probe provides for continuous monitoring of tissue blood flow (perfusion) in the target tissue or organ. Please see the Bowman Perfusion Monitor, Model 500, System Operating Manual for use instructions. Calibration data is stored within a memory chip in the QFlow 500 Probe connector. When the Probe is connected to the Bowman Perfusion Monitor, the system calibrates automatically, and is ready for use.  
This device should only be used by physicians with proper training in the area of application.  
Read and understand these instructions before beginning procedure.  
**CAUTION: Federal law (U.S.A.) restricts this device to sale or use by or on the order of a physician.**

**2. Indications**  
The QFlow 500 Probe is intended for extravascular monitoring of microcirculation blood flow in buried tissues. Examples of this application include but are not limited to: 1) the monitoring of buried muscle following three muscle transfer; 2) monitoring soft tissue microcirculation following reconstructive surgery such as facial reconstruction and; 3) monitoring cerebral blood flow during and following neurosurgery for head trauma.  
Use of product on other tissues has not been demonstrated. Demonstration is in progress.

**3. Contraindications**  
The QFlow 500 Probe is not intended for use other than described in **Indications**. Patients with contraindications for needle insertion into tissue, including those with coagulopathies and increased susceptibility to infection, should not be considered candidates for this device.

**4. Risks and complications**  
Maintenance of sterility during placement and subsequent handling is essential.  
Sterile technique should be used at all times when inserting, correcting or adjusting the Probe.

**5. Precautions**  
• Bending of the Probe 1 cm from the tip can result in damage and can impair performance.  
• Do not stretch or whip the Probe.  
• Exercise care when handling and inserting the Probe.  
• Use aseptic technique throughout the procedure.  
• Maintain the insertion site with regular redressing and use aseptic techniques.  
• This device takes a local (approximately 0.3 ml of tissue) measurement of blood flow, and regional flow extrapolation from the local blood flow measurement is not claimed.  
• This product may not be compatible with MRI diagnostic equipment; therefore, it is not recommended to conduct MRI imaging while this device is in use.

**6. Contents**  
One (1) sterile, QFlow 500 Probe, for single patient use.  
**CAUTION: Products labeled sterile are certified sterile unless package is opened or damaged.**

**7. Procedure**  
1. Check device package and contents for shipping or handling damage. Return any damaged or opened devices to Hemedex, Inc.  
2. Open sterile pouch and place device components into the sterile field.  
3. Connect the Probe to the Bowman Perfusion Monitor umbilical cord using aseptic technique.  
4. Establish an access track to the target tissue (or organ) with an introducer sheath, trocar or needle, and follow standard procedures for implanting Probe through the access track into the target tissue. Insert the Probe into the tissue such that at least 1 cm, up to the first colored band on the catheter, is buried within the tissue.  
5. See the System Operating Manual for the Bowman Perfusion Monitor, Model 500 for instructions on the use and operation of the monitor.  
**CAUTION: Do not leave the probe in place for more than 10 days.**

Parameter	Specification
Probe Diameter (O.D.)	1.15 mm (3.5 French or 18 gauge)
Perfusion Accuracy	±20 ml/min-100g
Probe Length	38 cm
Temperature Accuracy	±0.3 °C
Storage Time	Up to 10 days

Specifications are nominal values.  
**8. Cleaning, Disinfection, and Disposal**  
The QFlow 500 Probe is supplied sterile for single patient use only. Do not resterilize. If reused there is risk of infection or product failure. The single use Probe must be disposed of in an appropriate biohazard container in accordance with hospital practice.

**10. Warranty**  
This probe is warranted to be free from defects in material and workmanship for one year from the date of invoice and for a single patient use. Alterations or repairs done by persons not specifically authorized by Hemedex, Inc. will void this warranty. The foregoing warranty shall apply only to the original buyer. It is not valid for any anticipated profit, consequential damage or loss of income by the purchaser or use of this equipment.  
THE ABOVE WARRANTIES ARE IN LIEU OF ALL OTHER WARRANTIES, EITHER EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF FITNESS OR MERCHANTABILITY.

**11. Explanation of Symbols**

	Caution, consult accompanying documents		Manufacturer		Authorized EU Representative:
	Catalog Number		Manufacturer		Authorized EU Representative:
	Batch code		Manufacturer		Authorized EU Representative:
	Date of manufacture		Manufacturer		Authorized EU Representative:
	Use by		Manufacturer		Authorized EU Representative:

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