HEMEDIX
English
UMBILICAL CORD
For Use with the QFlow 500 Probe™ and Bowman Perfusion Monitor™
INSTRUCTIONS FOR USE

1. Introduction
The Umbilical Cord is used to connect the QFlow 500 Probe to the Bowman Perfusion Monitor, Model 500.
Please see the Bowman Perfusion Monitor, Model 500, System Operating Manual for use instructions.
This device should only be used by physicians with proper training in the area of application.
Read and understand these instructions before beginning procedure.
CAUTION: Federal law (U.S.A.) restricts this device to sale or use by or on the order of a physician.

2. Precautions
• This product may not be compatible with MRI diagnostic equipment; therefore, it is not recommended to conduct MRI imaging while this device is in use.
• Replace Umbilical Cord after 20 autoclave cycles or 200 days of use, whichever comes first.

3. Contents
One (1) non-sterile, reusable Umbilical Cord.

4. Procedure
1. Check device package and contents for shipping or handling damage. Return any damaged or opened devices to Hemedix, Inc.
2. Connect the Umbilical Cord to the QFlow 500 Probe using aseptic technique.
3. Connect the Umbilical Cord to the Bowman Perfusion Monitor, routing the cord to ensure that it does not present a hazard. Use endoased clips where necessary.
4. See the System Operating Manual for the Bowman Perfusion Monitor, Model 500 for instructions on the use and operation of the monitor.
5. Disconnect Umbilical Cord from Monitor by turning connector Locking Ring counterclockwise and then gently pulling Locking Ring away from Monitor.

5. Cleaning, Disinfection, and Disposal Umbilical Cord:
The umbilical cord must be cleaned and disinfected prior to each use. Blood and tissue may be removed from the cables and umbilical cord with a towel that has been dampened with an aqueous soap solution. Disinfectants containing a high percentage of alcohol or phenol may damage the cables.
If seal solutions, blood or other body fluids enter the cable connectors or sockets, the following cleaning procedure may be followed: The contaminated area of the connector may be rinsed with distilled water until all remains and contamination have been removed. A small brush can be used. The cord may be immersed for approximately 5 seconds in ethanol. The alcohol must be removed by the manufacturer's instructions.
Follow the autoclave chamber manufacturer's instructions. Do not overload the autoclave chamber as pockets may form that prevent steam penetration. Use the pre-vacuum cycle and autoclave as follows:
Temperature: 270 °F to 275 °F (132 °C to 135 °C)
Exposure Time: 3 to 4 minutes

IMPORTANT: It is recommended that the institution employ procedures that include the use of biological indicators in order to determine the effectiveness of the sterilization process.

6. Warranty
This umbilical cord is warranted to be free from defects in material and workmanship for one year from the date of invoice and for a single patient use. Alterations or repairs done by persons not specifically authorized by Hemedix, Inc. will void this warranty. The foregoing warranty shall apply only to the original buyer. It is not valid for any anticipated profits, consequential damages or a loss of time incurred by the buyer with the purchase or use of this equipment.

THE ABOVE WARRANTIES ARE IN LIEU OF ALL OTHER WARRANTIES, EITHER EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF FITNESS OR MERCHANTABILITY.

7. Explanation of Symbols
Caution, consult accompanying documents
Date of manufacture
Batch code
Manufacturer
Authorized European community representative
EC REP

Manufactured by: Hemedix, Inc. 222 Third Street, Suite 1123, Cambridge, MA 02142, USA. Tel: +31 (0)70 345 8570, Fax: +31 (0)70 346 7299.
Authorized EU Representative: EMERGO EUROPE, Molenaarsstraat 15, 2513 BH, The Hague, The Netherlands. Tel: +31 (0)70 345 8570, Fax: +31 (0)70 346 7299.
PIN H30990005, Rev. B

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